



Achilles Boot



Sizes: xs, s, m, l, xl

Our Hi Top boot with heel wedges places the foot in a plantar-flexed position using wedges in the heel. As the condition improves, wedges are removed, gently stretching the Achilles tendon and relieving tension on the tendon during rehab.

Common Examples of Use

- Immediate application following Achilles injury
- Reconstruction of Achilles tendon
- Sever's Disease

Product Features

- Mid-Calf or Standard height

Category: Walker Boots.

Part Numbers

Part Numbers

Part #	Description
AL03601XBB	Achilles Boot, S-XL
AL23601XBB	Achilles Boot, w/ Air, S-XL

Product Sizing Chart

Achilles Boot / Hi Top / Lo Top			
Size #	Size	Men's Shoe Size	Women's Shoe Size
X = 1	XS	1 – 2.5	1 – 3.5
X = 3	S	2.5 – 5.5	3.5 – 6
X = 5	M	5.5 – 9	6 – 9.5
X = 7	L	9 – 12.5	9.5 – 13
X = 9	XL	12.5 - 17	13+