



## Bunion Boot Walking Boot



Sizes: xs, s, m, l, xl

Pad-and-band holds the toe in position with constant gentle pressure. Soft Auto-mold sole cushions the foot while reducing unnecessary movement and flat innersole maintains neutral joint position.

### Common Examples of Use

- Post-operative of Hallux Abducto Valgus (HAV) surgery

### Product Features

- Pneumatic forefoot pad reduces swelling
- Bunion sock provided keeps great toe separate
- Fits left or right

Category: Walker Boots.

---

### Instructions for Use

[Instructions for Use](#)

### Part Numbers

## Part Numbers

Part #	Description
AL05800X	Bunion Boot, S-L
AL25800X	Bunion Boot w/ Air and Ankle Pad, S-L

## Product Sizing Chart

Bunion Boot Size Chart			
Size #	Size	Men	Women
X = 3	S	2.5 - 5.5	3.5 - 6
X = 5	M	5.5 - 9	6 - 9.5
X = 7	L	9 - 12.5	9.5 - 13