



Wee ROM Post-Op

A post-op knee brace with telescoping bars to accommodate varying leg lengths.



Common Examples of Use

- Non-operative or post-op procedures to the lower leg
- ACL, PCL, MCL and LCL injuries
- Meniscus
- Cartilage
- Juvenile Osteochondritis Dissecans
- Patella injuries

Product Features

- Universal sizing
- Telescoping bars adjust from 16" - 22" (41 - 56 cm)
- Full range-of-motion hinge adjustable from 10° - 110°
- Extends to malleolus

Categories: Post Op, Knee.

Part Numbers

Part Numbers

Part #	Description
EK090000	Wee ROM Post-Op, Universal
CE Marking	