

Suction pump FV series





Connectors and safety valve

Surgical suction pump model FV

The possibility to choose between a standalone flask of 1 litre and a maximum of 2 flask of 4 litres each, makes FUTURVAC one of the most versatile aspirators on market at when it comes to collecting the liquid waste, and therefore adaptable for all your needs.



Anti-static wheels and with brake

TECHNICAL DATA		
	FV series SM	FV series SH
Chassis :	Polyethylene moulded	
Engine pump :	Dry rotary piston vacuum pump	
Flow :	Flow 2.7 m ³ /h (45 l/min)	Flow 5.4 m ³ /h (90 l/min)
Max. Vacuum :	Vacuum regulator to 0-680 mmHg (91 kPa)	Vacuum regulator to 0-700 mmHg (93.5 kPa)
Noise Level :	58 dB	
Voltage :	230 V/50 hertz (under request)	
Power :	150 VA	350 VA
Weight :	22 Kg approximate	
Size :	53 X 43 X 80 cm	
Norms :	UNE EN 60601 / EN 61000 / UNE EN ISO 10079	
Models:	FV-1321SM : 2 flask of 3 litres + 1 of 1 litre FV-1303SM : 3 flask of 1 litre FV-1220SM : 2 flask of 3 litres FV-1202SM : 2 flask of 1 litre FV-1240SM : 2 flask of 4 litres	FV-1321SH : 2 flask of 3 litres + 1 of 1 litre FV-1303SH : 3 flask of 1 litre FV-1220SH : 2 flask of 3 litres FV-1202SH : 2 flask of 1 litre FV-1240SH : 2 flask of 4 litres

OPTIONALS SPARES

- FB-194A** Bacteria filter
- FV010018** Collecting flask of 3 litres (2 units)
- FV010019** Collecting flask of 1 litres (2 units)
- FV010300** Collecting flask of 3 litres, with top and connectors
- FV010100** Collecting flask of 3 litres, with top and connectors
- FV010400** Collecting flask of 4 litres, with top and connectors
- FV010052** Pedal ON/OFF with 1,5 meters of wire
- FV010068** Foot switch by air with tube of 2,15 m



Bactericide filter



Possibility of to incorporate an ON/OFF pedal



Tecnologia, 2, nau 7

08780, PALLEJA (Barcelona), SPAIN

Tel.: (+34) 93 334 01 12 - Fax: (+34) 93 440 25 64

www.ordisi.com - comercial@ordisi.com - export@ordisi.com

Reg. FVseries SP-EN/003

