

Warranty **3** years



## ALU ROLL / EASY ROLL

Hydraulic variable height :60 cm to 88 cm (without mattress)

OVERALL DIMENSIONS : Length = 211cm

Possible width = 50, 60, 70, 80 cm

Sleep platform 187 x 60cm

SWL : 200 kg

ASSISTED BACKREST by jack gas

Removable mattress, thickness 8 cm, tissue PVC M1



## EMERGENCY STETCHER 2052

Hydraulic variable height with bilateral control (56,5 to 89 cm out of mattress )

Out-of-everything dimensions : 210 x 90 cm – mattress 195 x 70 cm Memory foam

SWL: 230 kg

4 castors Ø 200 mm with central locking and 5th directionnal wheel

Radio-translucent stratified compact

Trendelenburg antitrendelenburg

maximum angle of 18°