

BED PACK Y-CARE TECH PLUS

Description



1 - EasyNursing: The H-frame has been specifically designed to be lighter, with no retention zone. With its full-width sleeping surface, patient handling and comfort are optimised. We must work with you to reduce MSD.

2 - EasyMov: The structure of the Y-Care healthcare bed, equipped with 150 mm castors and a compact and resistant low frame, guarantees the intensive mobility if needed



PACK - EQUIPMENT


Head & Foot boards	PPS Polypropylène
Side Rails	Half Side Rails Polypropylène

WEIGHT

Weight bare structure	110 kg
Total Weight - Pack	131 kg

SAFETY CHARGE

Maximum patient weight	195 kg
SWL	260 kg



PP-YC-TECH-PLUS-F : Detailed Technical Specifications

PACK - EQUIPMENT

Head & Foot boards	PPS Polypropylène
Side Rails	Half Side Rails Polypropylène

WEIGHT

Weight bare structure	110 kg
Total Weight - Pack	131 kg

DIMENSIONS

Length bare structure	216,5 cm
Width bare structure	102 cm
Length structure + H&F Boards + siderails	214,5 cm
Width structure + H&F Boards + siderails	102 cm
Length sleeping surface	200 cm
Width sleeping surface	90 cm
Variable height minimum	39 cm
Variable height maximum	87 cm
Standard bed extension length	20 cm

TYPE

Variable height type	Y-Type
Type of engine	Electric
Type of brake	Centralized brake

WEIGHING SYSTEM

Weighed (Yes/No)	No
------------------	----

MATERIAL & COATING

Tray Structure	White Epoxy Steel
Base Structure	Black Epoxy Steel
Frame	White Epoxy Steel
Variable height	Black Epoxy Steel
Bedbase / sleeping surface	Compact HPL

MATERIAL & COATING

Wheels	Thermoplastic rubber, grey non-marking
--------	--

BASE

Type of footing	Castors
Wheels diameter	150 mm

INCLINATIONS

Backrest	70 °
Legrest	15 °
Knee brake section	38 °
Knee brakes	14 °
Translation	11,5 cm
Trendelenburg	15 °
Reverse Trendelenburg	15 °

PROTECTION SIGN

Bed	IPX4
Engines	IP66

SAFETY CHARGE

Maximum patient weight	195 kg
SWL	260 kg

CONFIGURATION

Configuration	Siderails - H&F boards
Number of lifting pole support	2
Number of IV pole support	2

REGULATION

Medical device class	Class I
CE marking	CE
European regulation (UE 2017/745)	Yes

FIRE CLASSIFICATION

Motorization	V0
Compact HPL	M1
Laminated	M1
Melaminated panel	M3

OPTIONS

ELECTRIC KNEE BREAK WITH ELECTRIC FOOTREST (V20) FOR Y-CARE BED
Réf. OP-L-404-20-YC



PATIENT CONTROL SATELLITE REMOTE CONTROL FOR Y-CARE BED
Réf. OP-L-450-FPP



CAREGIVER SUPERVISOR PLACED AT THE FOOT OF THE BED WITH REMOTE CONTROL FOR Y-CARE BED
Réf. OP-L-466-SPVS-YC



BACKUP BATTERY FOR Y-CARE BED
Réf. OP-L-470-B-YC



NIGHTLIGHT PLACED UNDER THE SLEEPING PLATFORM FOR Y-CARE BED
Réf. OP-L-550-LU-YC



LOW POSITION INDICATOR
Réf. OP-L-550-P-BASSE



USB PLUG ON TECH AND TECH PLUS SIDE RAILS FOR PATIENT
Réf. OP-L-550-USB



ALARME FOR DISENGAGED BRAKES
Réf. OP-L-600-ALARME-FC



ELECTRIC OUTLET AT THE FOOT OF THE BED
Réf. OP-P-ELEC-YC



OPTIONS

BILATERAL CPR HANDLE OF BACKREST FOR Y-CARE BED
Réf. OP-L-314-2-YC



OPTIONS

DOUBLE CASTORS 150 MM (X4)
Réf. OP-L-601-150S-FCDG



CASTORS 150 MM INTEGRAL SMART (X4)
Réf. OP-L-601-150S-FCS



5TH CASTOR FOR Y-CARE BED
Réf. OP-L-650-5



OPTIONS

MOBILE HEAD-BED LOW
Réf. OP-L-700-TM



OPTIONS

HANDLE TO EASE GETTING OUT OF BED (THE PAIR)
Réf. OP-L-860-PSL-YC



MANDATORY ADDITIONAL SUPPORT TO "APP" SIDE RAILS FOR Y-CARE BEDS - SOLD WITHOUT APP
Réf. OP-SUP-APP-YC



OPTIONS

BEDDING SUPPORT FOR Y-CARE BED
Réf. OP-LX-SLI-YC



HOOKS FOR URINE BAG HOLDER ON RAIL (X8)
Réf. OP-LX-SPU-YC



OPTIONS

RADIO-TRANSLUSCENT BACKREST FOR Y-CARE BED
Réf. OP-L-300-RADIO



ACCESSORIES

FLEXIBLE SUPPORT FOR WIRED REMOTE CONTROL FOR Y-CARE BED
Réf. AC-LX-ST2H-YC



ACCESSORIES

SET OF CLINICA HEAD AND FOOT BOARDS - VERTICAL FITTING
Réf. AC-L-750-CLINICA-Y



ACCESSORIES

SET OF CLASSIC FOLDABLE SIDE RAILS - 3 BARS FOR BED Y4
Réf. AC-L-813-Y4



SET OF RAISED AND FOLDABLE SIDE RAILS 3/4 - 3 BARS
Réf. AC-L-813R-Y4



SET OF FOLDABLE STEEL SIDE RAILS 3/4
Réf. AC-L-814-Y4



SET OF RAISED AND FOLDABLE STEEL SIDE RAILS 3/4 - 4 BARS - ADAPTED TO GROWN CHILDREN
Réf. AC-L-814R-Y4



SET OF PADDED PROTECTIONS FOR CLASSIC SIDE RAILS (THE PAIR)
Réf. AC-LX-EVO07-2440P



ACCESSORIES

FIXED LIFTING POLE WITH STANDARD HANDLE
Réf. AC-LX-P1010



FIXED LIFTING POLE WITH AUTOMATIC REEL
Réf. AC-LX-P1010-ENR



LIFTING POLE SWIVELLED BY A PULL WITH STANDARD HANDLE - STEEL EPOXY
Réf. AC-LX-P1040



LIFTING POLE SWIVELLED BY A PULL WITH AUTOMATIC REEL
Réf. AC-LX-P1040-ENR



SUPPORT 4 PLASTIC HOOKS FOR LIFTING POLE
Réf. AC-LX-PCRO



ACCESSORIES

SET OF ADDITIONAL BARRIERS FOR CLASSIC SIDE RAILS (THE PAIR)
Réf. AC-LX-APP-Y4



ACCESSORIES

IV POLE WITH FIXED HEIGHT WITH 4 HOOKS
Réf. AC-LX-TPS-4C-Y4



IV POLE WITH HEIGHT ADJUSTABLE BY SCREW SYSTEM WITH 4 HOOKS
Réf. AC-LX-TPSR



IV POLE WITH VARIABLE HEIGHT BY SCREW SYSTEM - 4 HOOKS - STRAIGHT SHAPE
Réf. AC-LX-TPSR-2



ACCESSORIES

EXTENSION FOR WAFFLE FOAM MATTRESS + COVER - 20 X 90 X 14 CM
Réf. AC-LX-M315-Y4



MATTRESS + COVER - VISCOELASTIQUE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION
Réf. AC-LX-M325-Y4



ACCESSORIES

COVER FOR MATTRESS - WIDTH 90 CM - TRUNCATED ANGLES

Réf. AC-LX-H007-A-Y4



COVER FOR MATTRESS - WIDTH 90 CM - RIGHT ANGLES

Réf. AC-LX-H007-Y4



MATTRESS + COVER - COMFORT MULTI-ARTICULATION - SLEEPING 200 X 90 CM - TRUNCATED ANGLES

Réf. AC-LX-M207-A-Y4



MATTRESS + COVER - MULTI-ARTICULATION COMFORT - LYING SURFACE 200 X 90 CM - RIGHT ANGLES

Réf. AC-LX-M207-Y4



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION

Réf. AC-LX-M305-Y4



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 200 X 90 CM - TRUNCATED ANGLES

Réf. AC-LX-M407-A-Y4



MATTRESS + COVER - MULTI-JOINTS WAFFLE - LYING SURFACE 200 X 90 CM - RIGHT ANGLES

Réf. AC-LX-M407-Y4



MATTRESS + COVER - VISCOELASTIC - SLEEPING 200 X 90 CM - TRUNCATED ANGLES

Réf. AC-LX-M507-A-Y4



MATTRESS + COVER - VISCOELASTIC - LYING SURFACE 200 X 90 CM - RIGHT ANGLES

Réf. AC-LX-M507-Y4



MATTRESS + COVER - NON-MOTORIZED AIR - SLEEPING 200 X 90 CM

Réf. AC-LX-M707-Y4



MATTRESS + COVER - AUTOMATIC ALTERNATING LOW PRESSURE - SLEEPING 200 X 90 CM

Réf. AC-LX-M807



MATTRESS + COVER - PREVENT AND HELP THE TREATMENT OF PRESSURE SORES FOR VERY HIGH-RISK PATIENTS - SLEEPING 200 X 90 CM

Réf. AC-LX-M808



BED LOW EXTENSION FOR Y-CARE BED

Réf. AC-LX-ME-4-YC



ACCESSORIES

PATIENT PRESENCE AND/OR POSITION SENSOR IN THE BED

Réf. AC-LX-SANTINELLP



ACCESSORIES

STOPPER ON LOW BED FRAME OF Y-CARE BED - 2 POSITIONS

Réf. AC-LX-BU-YC



ACCESSORIES

URINAL HOLDER

Réf. AC-LX-ASUH



LABEL HOLDER

Réf. AC-LX-PET-YC



2 LITRES JAR HOLDER (INCOMPATIBLE WITH THE YMAGYN BED IN LOW POSITION)

Réf. AC-LX-SB2H



OXYGEN CYLINDER HOLDER FOR Y-CARE BED

Réf. AC-LX-SBAOH-YC



REDON HOLDER FOR 2 VIALS (INCOMPATIBLE WITH THE YMAGYN BED IN LOW POSITION)

Réf. AC-LX-SR2H



ACCESSORIES

MONITORING SUPPORT FOR Y-CARE BED

Réf. AC-LX-SMO-YC



SYRINGE PUMP SUPPORT FOR Y-CARE BED

Réf. AC-LX-SPS-YC



ACCESSORIES

TRACTION FRAME
Réf. AC-CADRE-TRAC-YC



KIT FOR CERVICAL SUSPENSION
Réf. AC-KIT-SUSP-CER-YC



KIT FOR PELVIC SUSPENSION
Réf. AC-KIT-SUSP-PEL-YC



TRACTION KIT FOR LOW LIMBS (KIT 1-A)
Réf. AC-KIT-TRAC-INF-YC



TRACTION KIT FOR UPPER LIMBS
Réf. AC-KIT-TRAC-SUP-YC

